Instructions for use –BEH\_ 19.07.2022 (Belotero® Hydro) RUSSIA

BELOTERO® HYDRO

MERZ AESTHETICS

## INSTRUCTIONS FOR USE FOR BELOTERO HYDRO

### Name of Medical Device

Belotero Hydro, sodium hyaluronate-based viscoelastic filler for rehydration of skin tissues, consisting of:

- implantable gel, in a 1.0 ml syringe - 1 pc.;
- label - 2 pcs.;
- Instructions for use – 1 pc.

### Intended Use

Belotero Hydro is designed only to be injected into the superficial dermis or the dermoepidermal junction for rehydration of skin.

### For use by medical practitioner.

Belotero Hydro should only be used by healthcare practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection. The treatment must be carried out under appropriate aseptic conditions.

### Description

Belotero Hydro is a sterile, non-pyrogenic, viscoelastic, colourless, transparent sodium hyaluronate gel of non-animal origin in a physiological phosphate citrate buffer containing glycerol.

### Presentation

Belotero Hydro is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe and two traceability labels.

### Composition

Sodium hyaluronate:	18 mg/ml
Glycerol:	21 mg/ml
Phosphate-citrate buffer pH 7 q.s:	gel volume

### Indications

Belotero Hydro is an injectable biodegradable implant indicated for filling of depressions due to facial atrophic scars and for skin rehydration and rejuvenation. Belotero Hydro enables to increase skin firmness and elasticity and to improve complexion and radiance of the face, neck, décolletage and hands.

### Posology and administration method

- Belotero Hydro is designed to be injected into the superficial dermis or the dermo-epidermal junction by a legally approved practitioner.
- Facial atrophic scar treatment is performed in the superficial part of the dermis in a single session.
- Skin rehydration and rejuvenation treatment is performed in the superficial dermis or dermo-epidermal junction with multipuncture or micro-papular injections in 3 sessions spaced 4 weeks apart. It is recommended to use the multipuncture technique over the whole face and perform a double series of parallel injections, the 2 series being perpendicular to each other. For the most dehydrated areas, it is recommended to use

## Instructions for use –BEH\_ 19.07.2022 (Belotero® Hydro) RUSSIA

the micro-papule technique consisting of repeated injections of small volumes. Patients should be informed that these micro-papules may be visible for 2 to 3 days.

- Inject Belotero Hydro slowly and apply the least amount of pressure necessary.
- In order to minimize the risks of potential complications, Belotero Hydro should only be used by healthcare practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection. The treatment must be carried out under appropriate aseptic conditions.
- Belotero Hydro must be injected into healthy, non-inflamed and previously appropriately disinfected skin.
- If necessary a local or loco-regional anaesthetic may be applied according to its instructions for use.
- It is recommended to use 30G1/2 needles with the product.
- Do not over-correct.
- There is no preclinical data available regarding the safety for injection of greater amounts of Belotero Hydro than 20 ml per year and 10 ml per session.
- Gently massage the treated area after the injection to distribute the product uniformly.
- Hyaluronic acid and glycerol are naturally occurring components of the skin and also responsible for water retention and balance in the tissues. Glycerol (polyalcohol) is a free scavenger which helps to stabilize sodium hyaluronate chains and the association of hyaluronic acid and glycerol should maintain the rehydrating effect up to 60 days, the product slowly degrading over this period.

### Contra-indications

Belotero Hydro is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate or to glycerol or to citric acid,
- In pregnant and breast-feeding women,
- In patients under 18 years old,
- In patients presenting a general infection.
- Do not inject Belotero Hydro into blood vessels. Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia or infarction.
- Do not inject Belotero Hydro into skin areas presenting active cutaneous inflammation or infection due to e.g. immunological, allergic, bacterial, fungal or viral causes.
- Do not inject Belotero Hydro into an area previously treated with a permanent dermal filler.
- Do not use Belotero Hydro in patient with joint, tendon, or vascular disease affecting the hand (when considering injection in the hand).

### Precautions for use

- Healthcare practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral haemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Practitioners should immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate healthcare practitioner specialist should an intravascular injection occur.

## Instructions for use –BEH\_ 19.07.2022 (Belotero® Hydro) RUSSIA

- In the absence of available clinical data on tolerance on the injection of Belotero Hydro in patients with previous medical history of or with an active autoimmune disease, or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject Belotero Hydro on a case-by-case basis depending on the nature of the disease as well as the associated treatment as it may worsen the existing patient health condition. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.
- No clinical data is available on the injection of Belotero Hydro into patient with a Fitzpatrick skin type VI.
- It is recommended not to inject Belotero Hydro in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.
- Belotero Hydro can be injected in the periorbital region only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area may be associated with an increased frequency and severity of side effects and adverse events.
- Belotero Hydro can be used in combination with other Belotero products during the same session only if injected in different facial areas. Instructions for use of each product should be followed.
- Belotero Hydro can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse®) only if injected in different facial areas. Practitioners should be experienced and patients appropriately selected as benefits but also adverse events can be cumulative and causality of adverse events could become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed. No clinical data are available on the injection of Belotero Hydro into an area already treated with other aesthetic products or procedures.
- Belotero Hydro must not be used in association with other aesthetic techniques such as peeling, or dermabrasion before complete healing of the last treatment. In any case, even if the healing occurs earlier, Belotero Hydro must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of Belotero Hydro with the above mentioned treatments.
- Belotero Hydro can be used in addition with laser treatments that should be conducted 1 or 2 months prior to hyaluronic acid injection. Patients using anti-coagulation, anti-platelet, or thrombolytic medications (e.g. warfarin), anti-inflammatory drugs (oral/injectable corticosteroids or non-steroidal anti-inflammatory drugs (NSAIDs; e.g. aspirin, ibuprofen), or other substances known to increase coagulation time (vitamins or herbal supplements, e.g. Vitamin E, garlic, ginkgo), from 10 days pre- to 3 days post injection may have increased reactions of hematomas, nodules or bleeding at the injection site.
- Injection of Belotero Hydro into persons with a history of previous herpetic eruption may be associated with herpes reactivation. Check the integrity of the inner packaging and the expiry date of the syringe prior to use. Do not use these products if the expiry date has lapsed or if the inner packaging has been opened or damaged.
- Do not transfer Belotero Hydro into another container and do not add other substances to the product.
- Only the gel is sterile, but not the outside of the syringe.
- Belotero Hydro must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an automatic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.
- Discard the syringe and the remaining product in the appropriate container after use.
- Do not re-sterilize and do not reuse due to the associated risks including infection.
- The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, peeling, Turkish baths and prolonged exposure to the sun or UV rays for at least 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

## Instructions for use –BEH\_ 19.07.2022 (Belotero® Hydro) RUSSIA

- If Belotero Hydro was injected with micro-papular technique, the patient can massage the papules only if the hands and the treated area have been previously disinfected.

### Warnings

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that Belotero Hydro does not come into contact with such substances.
- There is no known interaction with local or loco-regional anesthetics.

### Side effects and adverse events

Patients must be informed by the practitioner about possible side effects and adverse events before treatment.

- Side effects:

Injection site reactions may occur following injection into the skin but disappear spontaneously within a few days. This includes swelling, nodule or lump/bump, bruising, induration, erythema/redness, pain, discoloration and pruritus/itching. These injection site reactions are generally of mild or moderate intensity. A transient bleeding may also occur at the injection site and usually stops spontaneously as soon as the injection is finished.

- Adverse events:

In occasional cases, one or more of the following may occur either immediately or as a delayed reaction: inflammation, rash, burning sensation, pruritus/itching, hematoma, (lymph) edema, swelling, induration, lump/bump or nodule, granuloma, necrosis, vascular occlusion, hyper- or hypo-pigmentation, Tyndall effect, blister, and loss of performance. The risk of granuloma, ischemia or necrosis, and vascular occlusion is higher with deep injections and high volumes.

Rare cases of the following adverse events have been reported with hyaluronic acid products such as infection (e.g. cellulitis and abscess), scarring, hypersensitivity reaction or allergy to one of the product's components (e.g. hyaluronic acid, glycerol, citric acid), as well as sarcoid granuloma formation in subjects with hepatitis C and interferon treatment.

Isolated cases of visual impairment or blindness following unintentional intra-arterial injection have been reported in literature.

Patients should be instructed to report any side effects which last for more than one week and any adverse event as soon as it occurs to his/her practitioner, especially if patient has changes in his/her vision, signs of a stroke (including sudden difficulty speaking, numbness or weakness in his/her face, arms or legs, difficulty walking, face drooping, severe headache, dizziness, or confusion), white appearance of the skin, or unusual pain during or shortly after treatment. The practitioner may then refer the patient to the appropriate treatment.

Patients with lighter skin types are more likely to develop injection-related adverse events. However, patients with skin of color are more likely to develop post-inflammatory hyperpigmentation and/or hypertrophic scar/keloid formation following injection procedures. Patients with specific ethnic characteristics, e.g. Asian population, should be informed of a higher risk of tissue reactions, e.g. itching, swelling, erythema, inflammation.

### Storage

Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

## Instructions for use –BEH\_ 19.07.2022 (Belotero® Hydro) RUSSIA

### Disposal

Used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with medical procedures adopted in medical institutions and local, state or federal regulations.

### Requirements for environmental protection

The device should be disposed of in accordance with any applicable local/national regulations. With proper disposal, there is no harm to the environment.

### References

Updated documentation may be available from ANTEIS S.A., Switzerland.

Caution

Consult instructions for use.


Do not use if package is damaged.

Single use product. Do not re-use.

Open the blister by pulling the blister lid following the arrow.

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.

Temperature limit of storage: 2°C – 25°C.

Date of manufacture

Use-by date

Batch code

Manufacturer

Instructions for use –BEH\_ 19.07.2022 (Belotero® Hydro) RUSSIA

Manufacturer of Belotero Hydro:


**ANTEIS S.A.**  
18 Chemin des Aulx  
CH-1228 Plan-les-Ouates  
Geneva, Switzerland

**Distributed by:**


Authorized Representative in RF  
LLC «Merz Pharma»: 125124, Moscow, str. 3-ya Yamskogo polya, 2, bld. 26, 2nd floor,  
premises VII, room 29; tel. +7 (495) 653 8555; www.merz.ru

Date of the instructions for use: 2022-07-19

**A**

**B**

Backstop in the right position during injection

Seen exclusively for the legalization of the below signature of Mrs Barbara, Lise, Marine Humbert.


*The undersigned notary assumes no responsibility for the validity and content of the present document.*

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