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GB INSTRUCTIONS FOR USE LEAFLET FOR BELOTERO® BASIC

Description
 BELOTERO® Basic is a sterile, apyrogenic, viscoelastic, clear, colourless transparent, isotonic gel. BELOTERO® Basic consists of cross-linked sodium hyaluronate from bio fermentative origin.

Presentation
 BELOTERO® Basic is presented in a 1 ml glass syringe for single use, sterilised by moist heat by autoclave. Each syringe is supplied with 2 27G1/2 sterile needles (the product is in a combined packaging carrying the CE mark) and 2 traceability labels.

Composition
 Cross-linked sodium hyaluronate: 22,5 mg/ml
 Phosphate buffer pH 7 q.s.: 1 ml

Indications
 BELOTERO® Basic is an injectable implant indicated for use to increase the volume of skin tissue which enables:

- the correction of facial depressions, either due to injury or age-related (scarring, moderate or deep wrinkles, furrows),
- remodelling of the contours of the face by injection into facial areas such as the cheeks, chin, etc.
- increasing the volume of the lips (the border of the red part of the lip, or increasing roundness).

Dosage and method of administration
 BELOTERO® Basic is designed to be injected into the middle or deep dermis by a legally approved practitioner. For successful treatment it is essential that the practitioner has received prior training in the injection technique for increasing the volume of skin tissues.

BELOTERO® Basic must be injected into healthy, non-inflamed, disinfected skin. Screw the sterile 27G1/2 needle onto the Luer Lock tip, holding the Luer Lock. Do not connect other types of needle. Remove the protective needle sheath turning it in the same direction as the direction used to fix the needle. Start the injection pointing the needle cutting blade downwards.


BELOTERO® Basic may be injected into the middle or deep dermis using multi-puncture, linear retracing, blanket or fan distribution techniques. It is recommended that the linear retracing injection technique be used to fill depressions and to increase the volume of the lips. The needle is inserted over almost its entire length keeping closely parallel to the skin plane and then withdrawn slowly as the substance is deposited in the tissues. The blanket technique is preferred to change the contours of the face. This consists of a double series of parallel injections using the retracing technique, the two series of injections being administered perpendicular to each other in two different thicknesses of the skin tissue. If the needle becomes blocked and the injection pressure becomes too great, stop the injection and change the needle. Over-correction is not necessary. After the injection, massage gently in order to distribute the product uniformly. BELOTERO® Basic can be applied together with BELOTERO® Soft as BELOTERO® Soft should only be injected in the upper dermis after BELOTERO® Basic has been injected into the deep dermis.

Contra indications
 BELOTERO® Basic is contra indicated:

- in cases of hypersensitivity to hyaluronic acid,
- in cases of auto immune disease or diabetes,
- in cases of recurrent sore throat,
- in cases of acute joint rheumatism,
- in cases of Osler-rendu endocarditis,
- in patients pre-disposed to cheloid scarring,
- in pregnant or breast-feeding women,
- in young people less than 18 years old.

Do not inject BELOTERO® Basic into blood vessels. Do not inject BELOTERO® Basic into the periorbital region (rings around the eyes, eyelids, crow's feet, wrinkles) and do not use to correct superficial wrinkles. Do not use BELOTERO® Basic to increase the volume of the breasts. Do not use BELOTERO® Basic in the bones, tendons, ligaments or muscles.

BELOTERO® Basic must not be used in one treatment session in association with other aesthetic medicine techniques such as peeling, dermabrasion or any types of laser treatment, but at the earliest after complete healing of such treatment. In all cases, even if healing occurs earlier, BELOTERO® Basic must not be used earlier than 2 weeks after the last treatment. There is no experience for the combination of BELOTERO® Basic with the above mentioned treatments.

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Precautions for use

Before treatment the patient must be informed about the device, its contraindications and possible side effects.

Before use, confirm the status of the individual sterility protector and the expiry date both for the syringe and for the needle. Do not use these products if the expiry date has lapsed or if the sterility bag is opened or damaged.

Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO® Basic does not come into contact with such substances.

Sensitive skin may be pretreated with a local anaesthetic patch or cream. It should be noted that anaesthesia may cause local redness or hypersensitivity.

Do not transfer BELOTERO® Basic into another container or use an automated injection system to inject this product. Do not mix the product with other products.

BELOTERO® Basic is not to be injected subcutaneously.

Discard the syringe and remaining product after use. Do not resterilise and do not reuse.

No clinical data are available about injection of BELOTERO® Basic into an area which has already been treated with another filling product.

The patient must avoid applying make-up for at least 12 hours after treatment and avoid saunas and Turkish baths and prolonged exposure to sun or UV light for 2 weeks after the treatment. Patients should also avoid pressure and/or handling the area of the face which has been treated.

Side effects

Possible side effects exist and must be described to the patient before treatment. Very slight bleeding may occur during the injection although this disappears spontaneously as soon as the injection is finished. In occasional cases one or more of the following may occur either immediately or as a delayed reaction:

- local oedema or slight erythema sometimes associated with stinging in the treated area: these reactions may last for a week,
- allergy to hyaluronic acid,
- swelling or nodules at the injection site,
- discolouration of the treated area.

The patient should be instructed to report any side effect which lasts for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

Storage

Store between 2 and 25°C. Protect from light and freezing. Avoid shocks.

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Distribution by Merz Pharmaceuticals GmbH